



Update **_2_0_347A**

Clinical Reminders

**VA-BATTLEFIELD ACUPUNCTURE/BATTLEFIELD
ACUPRESSURE NOTE UPDATE**

Install Guide

January 2025

Product Development
Office of Information Technology
Department of Veterans Affairs

Table of Contents

INTRODUCTION	1
INSTALL DETAILS	2
INSTALL EXAMPLE	2
POST INSTALLATION	5

Introduction

Based on feedback from the field, the following changes have been made to the Battlefield Acupuncture Template

Here is a summary of the changes:

- 1) Addition of language for absolute contraindications in alignment with national training guidance
- 2) More detailed information in the note regarding consent and risks/benefits of the procedure(s)
- 3) Clarification of time statement to conform with coding guidance following consultation with HIMS
- 4) Update of aftercare instructions to match current training and Veteran facing aftercare handouts
- 5) The version is now 4.1

UPDATE_2_0_347A contains 1 Reminder Exchange entry:

UPDATE_2_0_347A VA-BATTLEFIELD ACUPUNCTURE/ACUPRESSURE
UPDATE

The exchange file contains the following components:

TIU TEMPLATE FIELD

VA-TEXT (3 CHAR)R
VA-BFA R-L-BL
WORD PROCESSING SHORT
BFA SUPPLEMENTAL QUESTIONS HYPERLINK

HEALTH FACTORS

VA-REMINDER UPDATES [C]
VA-UPDATE_2_0_347A
VA-HEALTH AND WELL BEING [C]
VA-WHS - BAA
VA-WHS - BFA
VA-BFA CATEGORY [C]
VA-BFA BAA TREATMENT INDIVIDUAL VISIT
VA-BFA BAA TREATMENT GROUP VISIT
VA-BFA ALL POINTS USED
VA-BFA NEEDLE TYPE OTHER
VA-BFA NEEDLE TYPE STANDARD NEEDLES
VA-BFA NEEDLE TYPE EAR SEED
VA-BFA NEEDLE TYPE EAR PRESS TACKS
VA-BFA NEEDLE TYPE SEMI PERM ASP
VA-BFA POST BFA PAIN RATING
VA-BFA PRE BFA PAIN RATING
VA-BFA BAA TYPE OF VISIT FOLLOW UP

VA-BFA BAA TYPE OF VISIT INITIAL
VA-BFA COMPLICATIONS
VA-BFA PAIN CONTRIBUTED TO STRESS PAST 24H RATING
VA-BFA PAIN INTERFERE MOOD PAST 24H RATING
VA-BFA PAIN INTERFERE SLEEP PAST 24H RATING
VA-BFA PAIN INTERFERE ACTIVITY PAST 24H RATING

REMINDER SPONSOR

OFFICE OF PATIENT CARE SERVICES

REMINDER TAXONOMY

VA-BAA CPT
VA-BFA CPT

REMINDER TERM

VA-REMINDER UPDATE_2_0_347A

REMINDER DIALOG

VA-BFA BATTLEFIELD ACUPUNCTURE NOTE

Install Details

This update is being distributed as a web host file. The address for the host file is:
https://REDACTED/UPDATE_2_0_347A.PRD

The file will be installed using Reminder Exchange, programmer access is not required.

Installation:

=====

This update can be loaded with users on the system. Installation will take less than 10 minutes.

Install Example

To Load the exchange file, use LWH. The URL is
https://REDACTED/UPDATE_2_0_347A.PRD

```
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE  Create Exchange File Entry  LHF  Load Host File
CHF  Create Host File             LMM  Load MailMan Message
CMM  Create MailMan Message       LR   List Reminder Definitions
DFE  Delete Exchange File Entry  LWH  Load Web Host File
IFE  Install Exchange File Entry  RI   Reminder Definition Inquiry
IH   Installation History         RP   Repack
Select Action: Next Screen// LWH  Load Web Host File
```

You should see a message that the file was loaded successfully.

Search and locate an entry titled UPDATE_2_0_347A VA-BATTLEFIELD ACUPUNCTURE/ACUPRESSURE UPDATE in reminder exchange.

+Item	Entry	Source	Date Packed
573	UPDATE_2_0_347A VA-BATTLEFIELD ACUPUNCTURE/ACU PRESSURE UPDATE		
574	UPDATE_2_0_348 VA-OSP CSRE UPDATE		
575	UPDATE_2_0_348 VA-OSP CSRE UPDATE		
576	UPDATE_2_0_349 VA-ANGIOEDEMA ORDER CHECK		
577	UPDATE_2_0_35 VA-BENEFICIARY		


```

+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE  Create Exchange File Entry      LHF  Load Host File
CHF  Create Host File                 LMM  Load MailMan Message
CMM  Create MailMan Message           LR   List Reminder Definitions
DFE  Delete Exchange File Entry       LWH  Load Web Host File
IFE  Install Exchange File Entry      RI   Reminder Definition Inquiry
IH   Installation History             RP   Repack
Select Action: Next Screen// IFE    Install Exchange File Entry
Enter a list or range of numbers (1-1009): 573

```

At the **Select Action** prompt, enter **IFE** for Install Exchange File Entry
 Enter the number that corresponds with your entry titled UPDATE_2_0_347A VA-BATTLEFIELD ACUPUNCTURE/BATTLEFIELD ACUPRESSURE NOTE UPDATE (in this example it is entry 573, it will vary by site) that is dated 10/08/2024.

```

Package Version: 2.0P65

Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
  VA-BFA BATTLEFIELD ACUPUNCTURE NOTE

REMINDER TERM
  VA-REMINDER UPDATE_2_0_347A

Keywords:

Components:

+      Enter ?? for more actions  >>>
IA   Install all Components      IS   Install Selected Component
Select Action: Next Screen// IA

```

At the **Select Action** prompt, type **IA** for **Install all Components** and hit enter.
 Any items that are NEW in any of the below categories, you should **INSTALL**
SKIP any TIU TEMPLATE FIELD if you are prompted to overwrite
OVERWRITE any HEALTH FACTORS if you are prompted to overwrite

You will then be prompted to install the Dialog component.

```
Packed reminder dialog: VA-BFA BATTLEFIELD ACUPUNCTURE NOTE [NATIONAL DIALOG]
Item  Seq.  Dialog Findings                                     Type  Exists
  1    1    VA-BFA BATTLEFIELD ACUPUNCTURE NOTE                   dialog  X
  2    5    VA-BFA BATTLEFIELD ACUPUNCTURE NOTE VERSION           element  X
        Finding: *NONE*
  3    10   VA-BFA GP INTIAL/FOLLOW UP                             group   X
        Finding: *NONE*
  4    10.5 VA-BFA GP INITIAL ONLY                               group   X
        Finding: VA-BFA BAA TYPE OF VISIT INITIAL (HEALTH FACTOR)
  5    10.10 VA-BFA GP FOLLOW UP ONLY                          group   X
        Finding: VA-BFA BAA TYPE OF VISIT FOLLOW UP (HEALTH FACTOR)
  6    PXR  PXR COMMENT                                           prompt  X
  7    15   VA-BLANK LINE DIALOG & NOTE                           element  X
        Finding: *NONE*
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage     QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// IA
```

At the **Select Action** prompt, type **IA** to install the dialog – **VA-BFA BATTLEFIELD ACUPUNCTURE NOTE**

Install reminder dialog and all components with no further changes: Y// ES....

```
Packed reminder dialog: VA-BFA BATTLEFIELD ACUPUNCTURE NOTE [NATIONAL DIALOG]
VA-BFA BATTLEFIELD ACUPUNCTURE NOTE (reminder dialog) installed from exchange fi
Item  Seq.  Dialog Findings                                     Type  Exists
  1    1    VA-BFA BATTLEFIELD ACUPUNCTURE NOTE                   dialog  X
  2    5    VA-BFA BATTLEFIELD ACUPUNCTURE NOTE VERSION           element  X
        Finding: *NONE*
  3    10   VA-BFA GP INTIAL/FOLLOW UP                             group   X
        Finding: *NONE*
  4    10.5 VA-BFA GP INITIAL ONLY                               group   X
        Finding: VA-BFA BAA TYPE OF VISIT INITIAL (HEALTH FACTOR)
  5    10.10 VA-BFA GP FOLLOW UP ONLY                          group   X
        Finding: VA-BFA BAA TYPE OF VISIT FOLLOW UP (HEALTH FACTOR)
  6    PXR  PXR COMMENT                                           prompt  X
  7    15   VA-BLANK LINE DIALOG & NOTE                           element  X
        Finding: *NONE*
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage     QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// Q
```

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

```
Package Version: 2.0P65

Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
  VA-BFA BATTLEFIELD ACUPUNCTURE NOTE

REMINDER TERM
  VA-REMINDER UPDATE_2_0_347A

Keywords:

Components:

+      + Next Screen  - Prev Screen  ?? More Actions  >>>
IA  Install all Components          IS  Install Selected Component
Select Action: Next Screen// Q
```

You will then be returned to this screen. At the **Select Action** prompt, type **Q**.

Install complete.

Post Installation

1. In CPRS or VistA, open/review the VA-BFA BATTLEFIELD ACUPUNCTURE NOTE template and confirm the version is 4.1.