



**Update\_2\_0\_430**

**Clinical Reminders**

**VA-MISC Update 27**

**Install Guide**

---

**April 2024**

Product Development  
Office of Information Technology  
Department of Veterans Affairs

**Table of Contents**

**INTRODUCTION ..... 1**

**PRE-INSTALLATION ..... 5**

**INSTALL DETAILS ..... 6**

**INSTALL EXAMPLE ..... 6**

**POST INSTALLATION ..... 27**

## **Introduction**

This update contains minor spelling corrections and other changes/edits to existing clinical reminder content, a full list of the changes can be found here:

### **REDACTED LINK**

**UPDATE\_2\_0\_430 contains 1 Reminder Exchange entry:**

UPDATE\_2\_0\_430 VA-MISC UPDATE 27

**The exchange file contains the following components:**

### **ROUTINES**

PXRMASU  
PXRMPDEM

### **MH TESTS AND SURVEYS**

FAST  
CDR  
GDS DEMENTIA

### **REMINDER GENERAL FINDINGS**

VIEW PROGRESS NOTE TEXT

### **TIU TEMPLATE FIELDS**

VA-WELLBEING SIGNS LINK  
VA-BLANK  
NUM REQ  
GEN WORD PROCESSING  
VA-SAH IA HEAD INJURY TEXT  
YES -OR- NO  
VA-SAH IA CURRENT MH TX  
YES/NO/NA  
TEXT (1-20 CHAR)  
TEXT (6 CHAR)  
TEXT (1-20 CHAR) REQ  
VA-NALOXONE PATIENT GUIDE  
VA-DATE ONLY R  
VA-TEXT (1-50 CHAR)R  
VA-ONC OFF REASON TXT  
DATE & TIME (REQ)  
VANOD PU DESCRIPTIONS  
VA-VAAES PREVENTION AND TREATMENT OF PI E-VERSION  
VA VANOD ACTIVITIES FOR PRESSURE INJURY IN CLC  
VA VANOD FLOW FOR SPINAL CORD  
VA VANOD FLOW FOR INPATIENT MENTAL HEALTH  
VA VANOD FLOW FOR HBPC PRESSURE INJURIES

VA VANOD FLOW FOR ACUTE PRESSURE INJURIES  
VA-WHOLE HEALTH LIBRARY

## HEALTH FACTORS

VA-M2VA CARE PLAN [C]  
VA-M2VA PAIN CARE PATHWAY MOTIVATIONAL INTERVENTION  
ADDICTION  
VA-M2VA PAIN CARE PATHWAY MOTIVATIONAL INTERVENTION PAIN  
CARE  
VA-REMINDER UPDATES [C]  
VA-UPDATE\_2\_0\_430  
REMINDER FACTORS [C]  
LIFE EXPECTANCY <6 MONTHS  
VA LIFE EXPECTANCY [C]  
VA LIMITED LIFE EXPECTANCY  
WH MAMMOGRAM [C]  
WH BILATERAL MASTECTOMY  
VA-WH BREAST CANCER RISK ASSESSMENT [C]  
VA-WH BCRA N/A MALE, NO ESTROGEN 5 YEARS  
VA-WH BCRA SERVICE ACT AREA NO  
LIFE EXPECTANCY < 1 YEAR  
LIFE EXPECTANCY < 6 MONTHS  
MENTAL HEALTH [C]  
SEVERE CHRONIC COGNITIVE IMPAIRMENT  
HOMELESSNESS SCREENING [C]  
VA-FOOD ALREADY RECEIVING ASSISTANCE  
REFERRED TO DIETICIAN FOOD  
REFERRED TO SOCIAL WORK FOOD  
NO REFERRALS REQUESTED FOOD  
VA-FOOD 12MO ACTUALLY RAN OUT OFTEN  
VA-IPVAP/SAH [C]  
VA-IPVAP/SAH LEG NOT JUSTICE INVOLVED  
VA-IPVAP/SAH LEG JUSTICE INVOLVED  
SKIN ASSESSMENT [C]  
VA-VAAES SKIN INSPECTION PU/PI NEW  
VAAES PI PERI WOUND ERYTHEMA NON-BLANCH  
PI PERI WOUND INTACT  
PI PERI WOUND USUAL FOR ETHNICITY  
PI PERI WOUND SKIN TEMP  
PI PERI WOUND FLUCTUANCE  
PI PERI WOUND INDURATION  
PI PERI WOUND CREPITUS  
PI PERI WOUND MACERATION  
PI PERI WOUND ERYTHEMA  
VA-VAAES PI STAGE 2 #1 WOUND BED TISSUE TYPE - HYPER-  
PIGMENTATION

VA-VAAES PI STAGE 2 #1 WOUND BED TISSUE TYPE - OTHER  
VA-VAAES PI STAGE 2 #1 WOUND BED TISSUE TYPE - RUPTURED BLISTER  
VA-VAAES PI STAGE 2 #1 WOUND BED TISSUE TYPE - INTACT BLISTER  
VA-VAAES PI STAGE 2 #1 WOUND BED TISSUE TYPE - RED  
VA-VAAES PI STAGE 2 #1 WOUND BED TISSUE TYPE - PINK  
VA-VAAES PI STAGE 2 #1 WOUND BED TISSUE TYPE - HYPO-PIGMENTATION  
VA-VAAES PI STAGE 2 #5 WOUND BED TISSUE TYPE - HYPER-PIGMENTATION  
VA-VAAES PI STAGE 2 #5 WOUND BED TISSUE TYPE - OTHER  
VA-VAAES PI STAGE 2 #5 WOUND BED TISSUE TYPE - RUPTURED BLISTER  
VA-VAAES PI STAGE 2 #5 WOUND BED TISSUE TYPE - INTACT BLISTER  
VA-VAAES PI STAGE 2 #5 WOUND BED TISSUE TYPE - RED  
VA-VAAES PI STAGE 2 #5 WOUND BED TISSUE TYPE - PINK  
VA-VAAES PI STAGE 2 #5 WOUND BED TISSUE TYPE - HYPO-PIGMENTATION  
SKIN HIGH RISK FACTORS [C]  
MED DEVICE PI STAGE 2 FIFTH  
DRAINAGE APPEARANCE VISCOUS  
DRAINAGE APPEARANCE SERO-PURULENT  
DRAINAGE APPEARANCE PURULENT  
DRAINAGE APPEARANCE SANGUINEOUS  
DRAINAGE APPEARANCE SEROSANGUINEOUS  
DRAINAGE APPEARANCE SEROUS  
VA-VAAES DRAINAGE COLOR RED  
VA-VAAES DRAINAGE COLOR AMBER  
DRAINAGE COLOR GREEN  
DRAINAGE COLOR TAN  
DRAINAGE COLOR COLORLESS  
DRAINAGE AMOUNT SCANT  
DRAINAGE AMOUNT HEAVY  
DRAINAGE AMOUNT MODERATE  
DRAINAGE AMOUNT SMALL  
WOUND DRAINAGE  
DRAINAGE AMOUNT NONE  
PRESSURE ULCER STAGE [C]  
STAGE II #5  
VA-VAAES PI STAGE 2 #4 WOUND BED TISSUE TYPE - HYPER-PIGMENTATION  
VA-VAAES PI STAGE 2 #4 WOUND BED TISSUE TYPE - OTHER  
VA-VAAES PI STAGE 2 #4 WOUND BED TISSUE TYPE - RUPTURED BLISTER  
VA-VAAES PI STAGE 2 #4 WOUND BED TISSUE TYPE - INTACT BLISTER  
VA-VAAES PI STAGE 2 #4 WOUND BED TISSUE TYPE - RED  
VA-VAAES PI STAGE 2 #4 WOUND BED TISSUE TYPE - PINK  
VA-VAAES PI STAGE 2 #4 WOUND BED TISSUE TYPE - HYPO-PIGMENTATION  
MED DEVICE PI STAGE 2 FOURTH  
STAGE II #4

VA-VAAES PI STAGE 2 #3 WOUND BED TISSUE TYPE - HYPER-PIGMENTATION  
VA-VAAES PI STAGE 2 #3 WOUND BED TISSUE TYPE - OTHER  
VA-VAAES PI STAGE 2 #3 WOUND BED TISSUE TYPE - RUPTURED BLISTER  
VA-VAAES PI STAGE 2 #3 WOUND BED TISSUE TYPE - INTACT BLISTER  
VA-VAAES PI STAGE 2 #3 WOUND BED TISSUE TYPE - RED  
VA-VAAES PI STAGE 2 #3 WOUND BED TISSUE TYPE - PINK  
VA-VAAES PI STAGE 2 #3 WOUND BED TISSUE TYPE - HYPO-PIGMENTATION  
MED DEVICE PI STAGE 2 THIRD  
STAGE II #3  
VA-VAAES PI STAGE 2 #2 WOUND BED TISSUE TYPE - HYPER-PIGMENTATION  
VA-VAAES PI STAGE 2 #2 WOUND BED TISSUE TYPE - OTHER  
VA-VAAES PI STAGE 2 #2 WOUND BED TISSUE TYPE - RUPTURED BLISTER  
VA-VAAES PI STAGE 2 #2 WOUND BED TISSUE TYPE - INTACT BLISTER  
VA-VAAES PI STAGE 2 #2 WOUND BED TISSUE TYPE - RED  
VA-VAAES PI STAGE 2 #2 WOUND BED TISSUE TYPE - PINK  
VA-VAAES PI STAGE 2 #2 WOUND BED TISSUE TYPE - HYPO-PIGMENTATION  
MED DEVICE PI STAGE 2 SECOND  
STAGE II #2  
MED DEVICE PI STAGE 2  
STAGE II  
VA-SUICIDE RISK ASSESSMENT [C]  
VA-SUICIDE ACTION IMPLEMENT APT TRACKING REMIND  
VA-NALOXONE [C]  
VA-NALOXONE PATIENT GUIDE  
VA-COMMON CARRIER NO ACCESS [C]  
VA-COMMON CARRIER NO ACCESS TO POV  
VA-COMMON CARRIER ACCESS TO POV  
LIFE EXPECTANCY < 5 YEARS  
VA-ONC PRCA CSPC M1 PATHWAY [C]  
VA-ONC PRCA CSPC M1 LV NO ABI OFF PATH RSN  
VA-VAAES LINE [C]  
VA-VAAES LINE ASSESS PERI #3 D/C  
VA-VAAES PERIPHERAL IV 3 D/C  
VA-VAAES LINE ASSESS PERI #3 INSERT  
VA-VAAES PERIPHERAL IV 3 INSERT  
VA-VAAES LINE ASSESS PERI #2 D/C  
VA-VAAES PERIPHERAL IV 2 D/C  
VA-VAAES LINE ASSESS PERI #2 INSERT  
VA-VAAES PERIPHERAL IV 2 INSERT  
VA-VAAES LINE ASSESS PERI #1 D/C  
VA-VAAES PERIPHERAL IV 1 D/C  
VA-VAAES LINE ASSESS PERI #1 INSERT  
VA-VAAES PERIPHERAL IV 1 INSERT

**REMINDER SPONSORS**

VA POST 9/11 TRANSITION AND CASE MANAGEMENT

NATIONAL CENTER FOR HEALTH PROMOTION AND DISEASE PREVENTION (NCP)  
MENTAL HEALTH AND BEHAVIORAL SCIENCE STRATEGIC GROUP  
MENTAL HEALTH SERVICES  
HOMELESSNESS OFFICE  
NATIONAL CLINICAL PRACTICE GUIDELINE COUNCIL  
OFFICE OF NURSING SERVICE  
OFFICE OF TELEHEALTH SERVICES STORE AND FORWARD TELEHEALTH

#### **REMINDER COMPUTED FINDINGS**

VA-ASU USER CLASS  
VA-IS INPATIENT

#### **REMINDER TAXONOMY**

VA-TERMINAL CANCER PATIENTS

#### **REMINDER TERMS**

VA-M2VA PAIN CARE PATHWAY MOTIVATIONAL INTERVENTION  
VA-REMINDER UPDATE\_2\_0\_430  
VA-WH BCRA NOT INDICATED  
VA-COGNITIVE IMPAIRMENT  
VA-USER CLASS FOR COGNITIVE IMPAIRMENT  
VA-VAAES IV LINE#3 DC  
VA-VAAES IV LINE#3  
VA-VAAES IV LINE#2 DC  
VA-VAAES IV LINE#2  
VA-VAAES IV LINE#1 DC  
VA-VAAES IV LINE#1

#### **REMINDER DEFINITIONS**

VA-BL NSG PIV #3 INSERT  
VA-BL NSG PIV #2 INSERT  
VA-BL NSG PIV #1 INSERT

#### **REMINDER DIALOGS**

VA-GP MISC UPDATE 27

## **Pre-Installation**

1. From the **Reminder Dialog Management** menu, Inquire/Print (INQ) and capture/save the following two (2) ELEMENTS :

ELEMENT: VAL-OI SOCIAL WORK REFERRAL

ELEMENT: VAL-OI DIETARY REFERRAL

2. Check the **Reminder Exchange** Install History (IH) and ensure your facility has installed the following (previously released) Updates:  
*(Post-install, you will re-map orders that you have mapped as an additional finding)*

- UPDATE\_2\_0\_229A VA-MILITARY SEXUAL TRAUMA SCREENING REMINDER
- UPDATE\_2\_0\_374 VA-COMMON CARRIER POV ACCESS
- UPDATE\_2\_0\_404 VA-WHS WELL-BEING SIGNS TEMPLATE UPDATE
- UPDATE\_2\_0\_417 VA-MISCELLANEOUS UPDATE 26
- UPDATE\_2\_0\_419 VA-MATERNITY CARE COORDINATION DIALOG UPDATE

## Install Details

This update is being distributed as a web host file. The address for the host file is:

**https://REDACTED/UPDATE\_2\_0\_430.PRD**

The file will be installed using Reminder Exchange, programmer access is not required.

Installation:

=====

This update can be loaded with users on the system. Installation will take less than 30 minutes.

## Install Example

To Load the Web Host File. Navigate to Reminder exchange in Vista

```

+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE  Create Exchange File Entry      LHF  Load Host File
CHF  Create Host File                 LMM  Load MailMan Message
CMM  Create MailMan Message           LR   List Reminder Definitions
DFE  Delete Exchange File Entry       LWH  Load Web Host File
IFE  Install Exchange File Entry      RI   Reminder Definition Inquiry
IH   Installation History             RP   Repack
Select Action: Next Screen// LWH  Load Web Host File ←
Input the URL for the .prd file: https://[REDACTED]/UPDATE_2_0_430.PRD

```

At the **Select Action:** prompt, enter **LWH** for Load Web Host File

At the **Input the url for the .prd file:** prompt, type the following web address:

**https://REDACTED/UPDATE\_2\_0\_430.PRD**

You should see a message at the top of your screen that the file successfully loaded.

```

https://[REDACTED]/UPDATE_2_0_430.PRD successfully loaded.

```

Search for and locate the reminder exchange entry titled UPDATE\_2\_0\_430 VA-MISC UPDATE 27.



+Item	Entry	Source	Date Packed
542	UPDATE_2_0_430 VA-MISC UPDATE 27		02/28/2024@11:19
543	UPDATE_2_0_44 VA-MHRRTP TEMPLATES		05/09/2018@14:06
544	UPDATE_2_0_45 VA-MATERNITY CARE COORDINATOR NOTE		02/08/2019@10:51
545	UPDATE_2_0_46 VA-MH TOBACCO USE SCREEN		08/14/2018@09:34
546	UPDATE_2_0_47 VA-HEPATITIS B IMMUNIZATIONS		08/15/2018@13:48
+    + Next Screen    - Prev Screen    ?? More Actions    >>>			
CFE	Create Exchange File Entry	LHF	Load Host File
CHF	Create Host File	LMM	Load MailMan Message
CMM	Create MailMan Message	LR	List Reminder Definitions
DFE	Delete Exchange File Entry	LWH	Load Web Host File
IFE	Install Exchange File Entry	RI	Reminder Definition Inquiry
IH	Installation History	RP	Repack
Select Action: Next Screen// IFE    Install Exchange File Entry			
Enter a list or range of numbers (1-732): 542			

At the **Select Action** prompt, enter **IFE** for Install Exchange File Entry  
Enter the number that corresponds with your entry titled **UPDATE\_2\_0\_430 VA-MISC UPDATE 27** (in this example it is entry 542, it will vary by site). The exchange file date is 02/28/2024@11:19.

Component	Category	Exists
Source: [REDACTED]		
Date Packed: 02/28/2024@11:19:33		
Package Version: 2.0P65		
Description:		
The following Clinical Reminder items were selected for packing:		
REMINDER DIALOG		
VA-GP MISC UPDATE 27		
REMINDER TERM		
VA-COGNITIVE IMPAIRMENT		
VA-WH BCRA NOT INDICATED		
VA-REMINDER UPDATE_2_0_430		
VA-M2VA PAIN CARE PATHWAY MOTIVATIONAL INTERVENTION		
REMINDER DEFINITION		
+    Enter ?? for more actions    >>>		
IA	Install all Components	IS    Install Selected Component
Select Action: Next Screen// <b>IS</b>		

**\*\*\*There is a new reminder sponsor associated with this install. You will need to do an IS to selectively install reminder sponsor VA POST 9/11 TRANSITION AND CASE MANAGEMENT. If you don't install this first, you will receive errors when installing the health factors. Please complete this next step prior to continuing! \*\*\***

+	Component	Category	Exists
	REMINDER SPONSOR		
151	VA POST 9/11 TRANSITION AND CASE MANAGEMENT		
152	NATIONAL CENTER FOR HEALTH PROMOTION AND DISEASE PREVENTION (NCP)		X
153	MENTAL HEALTH AND BEHAVIORAL SCIENCE STRATEGIC GROUP		X
154	MENTAL HEALTH SERVICES		X
155	HOMELESSNESS OFFICE		X
156	NATIONAL CLINICAL PRACTICE GUIDELINE COUNCIL		X
157	OFFICE OF NURSING SERVICE		X
158	OFFICE OF TELEHEALTH SERVICES STORE AND FORWARD TELEHEALTH		X
	REMINDER COMPUTED FINDINGS		
	VA-ASU USER CLASS		X
	VA-IS INPATIENT		X
+	+ Next Screen - Prev Screen ?? More Actions		>>>
IA	Install all Components	IS	Install Selected Component
Select Action: Next Screen// IS Install Selected Component			
Select Component(s): (1-174): 151			

At the **Select Action** prompt, type **IS** for Install Selected Component and hit enter.

At the Select Component, type 152 and hit enter

Select Action: Quit// **IS Install Selected Component**

Select Component(s): (1-175): **151**

REMINDER SPONSOR entry VA POST 9/11 TRANSITION AND CASE MANAGEMENT is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

+	Component	Category	Exists
	REMINDER SPONSOR		
152	VA POST 9/11 TRANSITION AND CASE MANAGEMENT		X
153	NATIONAL CENTER FOR HEALTH PROMOTION AND DISEASE PREVENTION (NCP)		X
154	MENTAL HEALTH AND BEHAVIORAL SCIENCE STRATEGIC GROUP		X
155	MENTAL HEALTH SERVICES		X
156	HOMELESSNESS OFFICE		X
157	NATIONAL CLINICAL PRACTICE GUIDELINE COUNCIL		X
158	OFFICE OF NURSING SERVICE		X
159	OFFICE OF TELEHEALTH SERVICES STORE AND FORWARD TELEHEALTH		X
	REMINDER COMPUTED FINDINGS		
	VA-ASU USER CLASS		X
	VA-IS INPATIENT		X
+    + Next Screen    - Prev Screen    ?? More Actions    >>>			
IA	Install all Components	IS	Install Selected Component
Select Action: Next Screen// <b>IA</b>			

At the **Select Action** prompt, type **IA** for Install all Components and hit enter.  
 Select Action: Next Screen// **IA Install all Components**

**You will see several prompts for components. The individual components for which you may (or may not) be prompted are listed below, along with the response you should provide.**

INSTALL	OVERWRITE	SKIP
<div style="border: 1px solid black; padding: 5px; width: fit-content;">           All new components         </div>	<div style="border: 1px solid black; padding: 5px;"> <p><b>TIU TEMPLATE FIELD:</b>            VANOD PU DESCRIPTIONS            VA-VAAE PREVENTTION AND TREATMENT OF PI E-VERSION            VA-VANOD ACTIVITIES FOR PRESSURE INJURY IN CLC            VA-VANOD FLOW FOR SPINAL CORD            VA-VANOD FLOW FOR INPATIENT MENTAL HEALTH            VA-VANOD FLOW FOR HBPC PRESSURE INJURIES            VA-VANOD FLOW FOR ACUTE PRESSURE INJURIES</p> <p><b>REMINDER TERM:</b>            VA-WH BCRA NOT INDICATED            VA-COGNITIVE IMPAIRMENT</p> <p><b>REMINDER DEFINITION:</b>            VA-BL NSG PIV #3 INSERT            VA-BL NSG PIV #2 INSERT            VA-BL NSG PIV #1 INSERT</p> </div>	<div style="border: 1px solid black; padding: 5px;"> <p><b>HEALTH FACTORS:</b>            LIFE EXPECTANCY &lt; 1 YEAR</p> <p><b>REMINDER TERM:</b>            VA-VAAE IV LINE#2 DC</p> </div>

*(The components outlined above are also listed out individually starting on the next page.)*  
**Any component NOT listed should be skipped.**

TIU TEMPLATE FIELD entry **VA-WELLBEING SIGNS LINK** is NEW,  
what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

TIU TEMPLATE FIELD entry named **VANOD PU DESCRIPTIONS** already exists  
but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// **Overwrite the current entry**

Are you sure you want to overwrite? N// YES

TIU TEMPLATE FIELD entry named **VA-VAES PREVENTION AND TREATMENT OF  
PI E-VERSION** already exists but the packed component is different,  
what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// **Overwrite the current entry**

Are you sure you want to overwrite? N// YES

TIU TEMPLATE FIELD entry named **VA VANOD ACTIVITIES FOR PRESSURE INJURY IN CLC** already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// **Overwrite the current entry**

Are you sure you want to overwrite? N// YES

TIU TEMPLATE FIELD entry named **VA VANOD FLOW FOR SPINAL CORD** already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// **Overwrite the current entry**

Are you sure you want to overwrite? N// YES

TIU TEMPLATE FIELD entry named **VA VANOD FLOW FOR INPATIENT MENTAL HEALTH** already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// **Overwrite the current entry**

Are you sure you want to overwrite? N// YES

TIU TEMPLATE FIELD entry named **VA VANOD FLOW FOR HBPC PRESSURE INJURIES** already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// **Overwrite the current entry**

Are you sure you want to overwrite? N// YES

TIU TEMPLATE FIELD entry named **VA VANOD FLOW FOR ACUTE PRESSURE INJURIES** already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// **Overwrite the current entry**

Are you sure you want to overwrite? N// YES

TIU TEMPLATE FIELD entry **VA-WHOLE HEALTH LIBRARY** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-M2VA CARE PLAN [C]** is NEW,  
what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-M2VA PAIN CARE PATHWAY MOTIVATIONAL INTERVENTION ADDICTION** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-M2VA PAIN CARE PATHWAY MOTIVATIONAL INTERVENTION PAIN CARE** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-UPDATE\_2\_0\_430** is NEW,  
what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**



HEALTH FACTORS entry named **LIFE EXPECTANCY < 1 YEAR** already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// **Skip, do not install this entry**

HEALTH FACTORS entry **VA-VAAES PI STAGE 2 #1 WOUND BED TISSUE TYPE -HYPER-PIGMENTATION** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-VAAES PI STAGE 2 #1 WOUND BED TISSUE TYPE - OTHER** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-VAAES PI STAGE 2 #1 WOUND BED TISSUE TYPE - RUPTURED BLISTER** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**



HEALTH FACTORS entry VA-VAAES PI STAGE 2 #1 WOUND BED TISSUE TYPE - INTACT BLISTER is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// Install

HEALTH FACTORS entry VA-VAAES PI STAGE 2 #1 WOUND BED TISSUE TYPE - RED is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// Install

HEALTH FACTORS entry VA-VAAES PI STAGE 2 #1 WOUND BED TISSUE TYPE - PINK is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// Install

HEALTH FACTORS entry VA-VAAES PI STAGE 2 #1 WOUND BED TISSUE TYPE - HYPO-PIGMENTATION is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// Install

HEALTH FACTORS entry **VA-VAES PI STAGE 2 #5 WOUND BED TISSUE TYPE -HYPER-PIGMENTATION** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-VAES PI STAGE 2 #5 WOUND BED TISSUE TYPE - OTHER** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-VAES PI STAGE 2 #5 WOUND BED TISSUE TYPE - RUPTURED BLISTER** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-VAES PI STAGE 2 #5 WOUND BED TISSUE TYPE - INTACT BLISTER** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-VAES PI STAGE 2 #5 WOUND BED TISSUE TYPE - RED** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-VAES PI STAGE 2 #5 WOUND BED TISSUE TYPE - PINK** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-VAES PI STAGE 2 #5 WOUND BED TISSUE TYPE - HYPO-PIGMENTATION** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-VAES PI STAGE 2 #4 WOUND BED TISSUE TYPE -HYPER-PIGMENTATION** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry VA-VAAES PI STAGE 2 #4 WOUND BED TISSUE TYPE - OTHER is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// Install

HEALTH FACTORS entry VA-VAAES PI STAGE 2 #4 WOUND BED TISSUE TYPE - RUPTURED BLISTER is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// Install

HEALTH FACTORS entry VA-VAAES PI STAGE 2 #4 WOUND BED TISSUE TYPE - INTACT BLISTER is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// Install

HEALTH FACTORS entry VA-VAAES PI STAGE 2 #4 WOUND BED TISSUE TYPE - RED is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// Install

HEALTH FACTORS entry **VA-VAES PI STAGE 2 #4 WOUND BED TISSUE TYPE - PINK** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-VAES PI STAGE 2 #4 WOUND BED TISSUE TYPE - HYPO-PIGMENTATION** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-VAES PI STAGE 2 #3 WOUND BED TISSUE TYPE -HYPER-PIGMENTATION** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-VAES PI STAGE 2 #3 WOUND BED TISSUE TYPE - OTHER** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry VA-VAAES PI STAGE 2 #3 WOUND BED TISSUE TYPE - RUPTURED  
BLISTER is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// Install

HEALTH FACTORS entry VA-VAAES PI STAGE 2 #3 WOUND BED TISSUE TYPE - INTACT  
BLISTER is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// Install

HEALTH FACTORS entry VA-VAAES PI STAGE 2 #3 WOUND BED TISSUE TYPE - RED is  
NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// Install

HEALTH FACTORS entry VA-VAAES PI STAGE 2 #3 WOUND BED TISSUE TYPE - PINK is  
NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// Install

HEALTH FACTORS entry **VA-VAAES PI STAGE 2 #3 WOUND BED TISSUE TYPE - HYPO-PIGMENTATION** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-VAAES PI STAGE 2 #2 WOUND BED TISSUE TYPE -HYPER-PIGMENTATION** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-VAAES PI STAGE 2 #2 WOUND BED TISSUE TYPE - OTHER** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-VAAES PI STAGE 2 #2 WOUND BED TISSUE TYPE - RUPTURED BLISTER** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-VAAES PI STAGE 2 #2 WOUND BED TISSUE TYPE - INTACT BLISTER** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-VAAES PI STAGE 2 #2 WOUND BED TISSUE TYPE - RED** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-VAAES PI STAGE 2 #2 WOUND BED TISSUE TYPE - PINK** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-VAAES PI STAGE 2 #2 WOUND BED TISSUE TYPE - HYPO-PIGMENTATION** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**



HEALTH FACTORS entry **VA-SUICIDE ACTION IMPLEMENT APT TRACKING REMIND** is NEW,  
what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

HEALTH FACTORS entry **VA-NALOXONE PATIENT GUIDE** is NEW,  
what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

REMINDER TERM entry **VA-M2VA PAIN CARE PATHWAY MOTIVATIONAL INTERVENTION** is  
NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

REMINDER TERM entry **VA-REMINDER UPDATE\_2\_0\_430** is NEW,  
what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// **Install**

REMINDER TERM entry named VA-WH BCRA NOT INDICATED already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- M Merge findings
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// Overwrite the current entry

Are you sure you want to overwrite? N// YES

REMINDER TERM entry named VA-COGNITIVE IMPAIRMENT already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- M Merge findings
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// Overwrite the current entry

Are you sure you want to overwrite? N// YES

REMINDER TERM entry named VA-VAES IV LINE#2 DC already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- M Merge findings
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// Skip, do not install this entry

REMINDER DEFINITION entry named VA-BL NSG PIV #3 INSERT already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// Overwrite the current entry

Are you sure you want to overwrite? N// YES

REMINDER DEFINITION entry named VA-BL NSG PIV #2 INSERT already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// Overwrite the current entry

Are you sure you want to overwrite? N// YES

REMINDER DEFINITION entry named VA-BL NSG PIV #1 INSERT already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// Overwrite the current entry

Are you sure you want to overwrite? N// YES

**You will be prompted to install the dialog component:**

```
Packed reminder dialog: VA-GP MISC UPDATE 27
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-GP MISC UPDATE 27 Finding: *NONE*	group	
2	10	VA-MISC 27 TEMPLATE FIELDS Finding: *NONE*	element	
3	20	VA-TXT TRI IOP DONE-DATE Finding: *NONE*	element	X
4	30	VA-ONC PRCA CSPC M1 LOW CAND ABI NO OFF PATH REASON Finding: VA-ONC PRCA CSPC M1 LV NO ABI OFF PATH RSN (HEALTH FACTOR)	element	X
5		VA-ONC OFF PATHWAY	prompt	X

```
+      + Next Screen  - Prev Screen  ?? More Actions
DD Dialog Details    DT Dialog Text    IS Install Selected
DF Dialog Findings   DU Dialog Usage   QU Quit
DS Dialog Summary    IA Install All
Select Action: Next Screen// IA
```

At the **Select Action** prompt, type **IA** to install the dialog – **VA-GP MISC UPDATE 27**

Select Action: Next Screen// **IA Install All**

Install reminder dialog and all components with no further changes: **Y// ES**

```
Packed reminder dialog: VA-GP MISC UPDATE 27
VA-GP MISC UPDATE 27 (group) installed from exchange file.
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-GP MISC UPDATE 27 Finding: *NONE*	group	X
2	10	VA-MISC 27 TEMPLATE FIELDS Finding: *NONE*	element	X
3	20	VA-TXT TRI IOP DONE-DATE Finding: *NONE*	element	X
4	30	VA-ONC PRCA CSPC M1 LOW CAND ABI NO OFF PATH REASON Finding: VA-ONC PRCA CSPC M1 LV NO ABI OFF PATH RSN (HEALTH FACTOR)	element	X
5		VA-ONC OFF PATHWAY	prompt	X

```
+      + Next Screen  - Prev Screen  ?? More Actions
DD Dialog Details    DT Dialog Text    IS Install Selected
DF Dialog Findings   DU Dialog Usage   QU Quit
DS Dialog Summary    IA Install All
Select Action: Next Screen// Q
```

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

Component	Category	Exists
Source: [REDACTED]		
Date Packed: 02/28/2024@11:19:33		
Package Version: 2.0P65		
Description:		
The following Clinical Reminder items were selected for packing:		
REMINDER DIALOG		
VA-GP MISC UPDATE 27		
REMINDER TERM		
VA-COGNITIVE IMPAIRMENT		
VA-WH BCRA NOT INDICATED		
VA-REMINDER UPDATE 2_0_430		
VA-M2VA PAIN CARE PATHWAY MOTIVATIONAL INTERVENTION		
REMINDER DEFINITION		
+ Enter ?? for more actions >>>		
IA	Install all Components	IS Install Selected Component
Select Action: Next Screen// Q		

You will be returned to this screen. Type Q to quit.

Install Completed

## Post Installation

1. Refer to the Inquire/Print (INQ) you saved during the Pre-Installation steps for the following two (2) ELEMENTS:

*(Map the quick orders identified in the pre-installation step to the element's additional finding)*

Element: **VAL-OI SOCIAL WORK REFERRAL**

FINDING ITEM:

Additional findings:

Select ADDITIONAL FINDING: Enter "Q.", then name of local consult quick order

Element: **VAL-OI DIETARY REFERRAL**

FINDING ITEM:

Additional findings:

Select ADDITIONAL FINDING: Enter "Q.", then name of local consult quick order

2. Locate the group: **VAL-GP MCC ALCOHOL USE SCREEN CR MAIN**

- a. In VISTA, navigate to the Reminder Dialog Management. Select Reminder Dialogs (DI), then Change View (CV) to Dialog Group View
- b. Locate the group titled: **VAL-GP MCC ALCOHOL USE SCREEN CR MAIN**
- c. Replace the group in Sequence 10 to be: **VA-ALCOHOL SCREENING VERSION**
- d. Remove/Delete Sequence 20: **VAL-GP ALC SCREEN AUDC**  
*(note: leave the other sequences as-is)*

```

NAME: VAL-GP MCC ALCOHOL USE SCREEN CR MAIN
DISABLE: <enter>
CLASS: LOCAL// <enter>
SPONSOR: <enter>
REVIEW DATE: <enter>
RESOLUTION TYPE: <enter>
ORDERABLE ITEM: <enter>
FINDING ITEM: <enter>
Additional findings: <enter>none
Select ADDITIONAL FINDING: <enter>
GROUP CAPTION: <enter>
PUT A BOX AROUND THE GROUP: <enter>
SHARE COMMON PROMPTS: <enter>
MULTIPLE SELECTION: <enter>
HIDE/SHOW GROUP: <enter>
GROUP HEADER DIALOG TEXT:
Edit? NO// <enter>
GROUP HEADER ALTERNATE P/N TEXT:
Edit? NO// <enter>
EXCLUDE FROM PROGRESS NOTE: YES// <enter>
SUPPRESS CHECKBOX: <enter>
NUMBER OF INDENTS: <enter>
INDENT PROGRESS NOTE TEXT: <enter>
Select SEQUENCE: 10
DIALOG ELEMENT: VA-TEXT ALCOHOL MISUSE SCREEN INSTRUCTIONS << Replace
// VA-ALCOHOL SCREENING VERSION << with
...OK? Yes//
EXCLUDE FROM PN TEXT: <enter>
Select SEQUENCE: 20 VAL-GP ALC SCREEN AUDC
SEQUENCE: 20// @
SURE YOU WANT TO DELETE THE ENTIRE '20' SEQUENCE? Y
Select SEQUENCE: <enter>
Select BRANCHING LOGIC SEQUENCE: <enter>
Checking reminder dialog for errors
NO ERRORS FOUND
Input your edit comments.
Edit? NO//<enter>

```

The group VAL-GP MCC ALCOHOL USE SCREEN CR MAIN will now have the following sequences when you complete this step:

- 10 VA-ALCOHOL SCREENING VERSION
- 30 TEXT BLANK LINE WITH TEMPLATE FIELD
- 40 VA-TEXT SUD CPG URL
- 50 TEXT BLANK LINE WITH TEMPLATE FIELD