



Update_2_0_460

Clinical Reminders

VA-VAAES SKIN INSPECTION/ASSESSMENT UPDATE

Install Guide

August 2024

Product Development
Office of Information Technology
Department of Veterans Affairs

Table of Contents

INTRODUCTION	1
PRE-INSTALLATION.....	14
INSTALL DETAILS.....	14
INSTALL EXAMPLE	14
POST INSTALLATION	18

Introduction

The Office of Nursing Informatics (ONI) strives to build national templates centered around reduced documentation burden, increased evidence-based awareness at point of care, and to support the ability to facilitate the creation of near real time reports on nursing sensitive indicators. The changes to this template were backend nursing informatics only:

Updates include:

Incorporating additional “group boxes” into template design (see screenshots below)

Update version number to 4.1

Added additional health factors to items that did not previously contain health factors.

There are no changes to content for the end users. The content aligns with the latest evidenced based practice guidelines. An update of nationally published guidelines aimed at preventing pressure ulcers/injury was developed by this VA workgroup resulting in a VHA Handbook 1180.02 in 2011 to a directive (VHA Directive 1352 Prevention and Management of Pressure Injuries, 3/21/ 2019). These quality indicators provide valuable data to further support the best possible patient care outcomes. The purpose of the national template is to preserve the ability to collect pertinent data for quality-driven healthcare. The content will not be altered. However, additional content may be added at local discretion. A yearly review will be conducted by the Office of Nursing Informatics to optimize core content and bring new evidence-based knowledge to our clinicians at point of care.

UPDATE_2_0_460 contains 1 Reminder Exchange entry:

UPDATE_2_0_460 VA-VAAES SKIN INSPECTION/ASSESSMENT UPDATE

The exchange file contains the following components:

REMINDER GENERAL FINDINGS

VIEW PROGRESS NOTE TEXT

REMINDER DIALOG LINK TYPE

VA-VANOD HF PU PROTOCOL NA NO RISK

TIU TEMPLATE FIELD

COMMENT

ED BOX

ED BLANK 3

ED BLANK 5

DATE/TIME1

VA-VANOD BLANK

TEXT (1-20 CHAR)
TEXT (6 CHAR)
VA-NSG FLAP PULSE
VA-NSG FLAP COLOR
VA-NSG FLAP TEXT
VA-NSG FLAP CAP
VA-WORD 1 LINE
VA-NSG FLAP TEMP
WORD 1 LINE
TEXT (1-20 CHAR) REQ
TEXT (1-30 CHAR) REQ
VANOD PU DESCRIPTIONS
WORD PROCESS (REQUIRED)
VA-WORD
TEXT (1-60 CHAR)
VA VANOD PRESSURE INJURY RISK FACTORS
BLANK TEXT FIELD FOR DIALOGS
VA VANOD ONS SHAREPOINT LINK
VA-VAAES PREVENTION AND TREATMENT OF PI E-VERSION
VA-VAAES PREVENTION AND TREATMENT OF PI
VA VANOD ACTIVITIES FOR PRESSURE INJURY IN CLC
VA VANOD FLOW FOR SPINAL CORD
VA VANOD FLOW FOR INPATIENT MENTAL HEALTH
VA VANOD FLOW FOR HBPC PRESSURE INJURIES
VA VANOD FLOW FOR ACUTE PRESSURE INJURIES
VA VANOD VHA DIRECTIVE 1352

HEALTH FACTORS

VA-REMINDER UPDATES [C]
VA-UPDATE_2_0_460
VA-VAAES URO [C]
VA-VAAES URO ADD OST PT PERF OTHER
VA-VAAES URO ADD OST PT PERF TOT ASST
VA-VAAES URO ADD OST PT PERF ASSIST
VA-VAAES URO ADD OST PT PERF IND
VA-VAAES URO OSTOMY MGN HYGIENE COMP
VA-VAAES URO OSTOMY MGN WAFER CHANGED
VA-VAAES URO OSTOMY MGN POUCH CHANGED
VA-VAAES GASTRO [C]
VA-VAAES GASTRO SUP OST PT PERF ASST DEV
VA-VAAES GASTRO SUP OST PT PERF TOT ASST
VA-VAAES GASTRO SUP OST PT PERF ASSIST
VA-VAAES GASTRO SUP OST PT PERF IND
VA-VAAES GASTRO OSTOMY MGN HYGIENE COMP

VA-VAAES GASTRO OSTOMY MGN WAFER CHANGED
VA-VAAES GASTRO OSTOMY MGN POUCH CHANGED
SKIN ASSESSMENT [C]
VA-VAAES SKIN/WOUND DRSG OPEN TO AIR
VA-VAAES SKIN/WOUND DRSG NEG PRESSURE
VA-VAAES SKIN/WOUND DRSG CHG PER ORDER
VA-VAAES SKIN/WOUND DRSG REINFORCED
VA-VAAES SKIN/WOUND DRSG INTACT
VA-VAAES SKIN/WOUND DRSG DRY
VA-VAAES SKIN/WOUND DRSG CLEAN
VA-VAAES SKIN/WOUND DRSG LOCATION 5
VA-VAAES SKIN/WOUND DRSG LOCATION 4
VA-VAAES SKIN/WOUND DRSG LOCATION 3
VA-VAAES SKIN/WOUND DRSG LOCATION 2
VA-VAAES SKIN/WOUND DRSG LOCATION 1
VANOD SKIN REMINDER TERMS [C]
VANOD UNABLE TO ASSESS OTHER
PATIENT REFUSAL END OF LIFE
PATIENT REFUSAL PERSONAL
PATIENT BEHAVIORALLY UNSTABLE
PATIENT MEDICALLY UNSTABLE
PATIENT OFF FLOOR
VANOD UNABLE TO ASSESS
VA-VAAES SKIN TELEHEALTH/VIRTUAL
SKIN TEMPERATURE
SKIN TURGOR
SKIN MOISTURE
SKIN COLOR
VA-VAAES URO ADD OST OTHER
VA-VAAES URO ADD OST ILEAL CONDUIT
VA-VAAES URO ADD OST CYSTOSTOMY
VA-VAAES URO ADD OST UROSTOMY
VA-VAAES URO OST LOCATION
VA-VAAES URO OST STOMA EST
VA-VAAES URO OST STOMA NEW
VA-VAAES URO ADD OST
VA-VAAES GASTRO SUP OST PT PERF OTHER
VA-VAAES GASTRO SUP OST OTHER
VA-VAAES GASTRO SUP OST ILEOSTOMY
VA-VAAES GASTRO SUP OST COLOSTOMY
VA-VAAES GASTRO OST LOCATION
VA-VAAES GASTRO SUP OST STOMA EST
VA-VAAES GASTRO SUP OST STOMA NEW
VA-VAAES GASTRO SUP OST OSTOMY

VANOD SKIN REASSESSMENT
VANOD SKIN INITIAL
VA-VAAES SKIN/WOUND DRSG NON-ADHERING
VA-VAAES SKIN/WOUND DRSG TRANSPARENT
VA-VAAES SKIN/WOUND DRSG FOAM DRSG
VA-VAAES SKIN/WOUND DRSG GAUZE WRAP
VA-VAAES SKIN/WOUND DRSG GAUZE
VA-VAAES SKIN/WOUND DRSG ANTI-MC
VA-VAAES INCISION/FLAP [C]
VA-VAAES INCISION/FLAP DRAIN VISCOUS
VA-VAAES INCISION/FLAP DRAIN SERO-PURUL
VA-VAAES INCISION/FLAP DRAIN PURULENT
VA-VAAES INCISION/FLAP DRAIN SANG
VA-VAAES INCISION/FLAP DRAIN SEROSANG
VA-VAAES INCISION/FLAP DRAIN SEROUS
VA-VAAES INCISION/FLAP DRAIN RED
VA-VAAES INCISION/FLAP DRAIN AMBER
VA-VAAES INCISION/FLAP DRAIN GREEN
VA-VAAES INCISION/FLAP DRAIN TAN
VA-VAAES INCISION/FLAP DRAIN COLORLESS
VA-VAAES INCISION/FLAP DRAIN SCANT
VA-VAAES INCISION/FLAP DRAIN LARGE
VA-VAAES INCISION/FLAP DRAIN MODERATE
VA-VAAES INCISION/FLAP DRAIN SMALL
VA-VAAES INCISION/FLAP DRAINAGE YES
VA-VAAES INCISION/FLAP DRAINAGE NO
VA-VAAES INCISION/FLAP #4 TEMP CHANGE
VA-VAAES INCISION/FLAP #4 MACERATION
VA-VAAES INCISION/FLAP #4 INDURATION
VA-VAAES INCISION/FLAP #4 ERYTHEMA
VA-VAAES INCISION/FLAP #4 INTACT
VA-VAAES INCISION/FLAP #4 OPEN/DEHISCED
VA-VAAES INCISION/FLAP #4 UNABLE
VA-VAAES INCISION/FLAP FLAP #4 IN PLACE
VA-VAAES INCISION/FLAP #4 WELL APPROX
VA-VAAES INCISION/FLAP #4 ASSESSED
VA-VAAES INCISION/FLAP #3 TEMP CHANGE
VA-VAAES INCISION/FLAP #3 MACERATION
VA-VAAES INCISION/FLAP #3 INDURATION
VA-VAAES INCISION/FLAP #3 ERYTHEMA
VA-VAAES INCISION/FLAP #3 INTACT
VA-VAAES INCISION/FLAP #3 OPEN/DEHISCED
VA-VAAES INCISION/FLAP #3 UNABLE
VA-VAAES INCISION/FLAP FLAP #3 IN PLACE

VA-VAAES INCISION/FLAP #3 WELL APPROX
VA-VAAES INCISION/FLAP #3 ASSESSED
VA-VAAES INCISION/FLAP #2 TEMP CHANGE
VA-VAAES INCISION/FLAP #2 MACERATION
VA-VAAES INCISION/FLAP #2 INDURATION
VA-VAAES INCISION/FLAP #2 ERYTHEMA
VA-VAAES INCISION/FLAP #2 INTACT
VA-VAAES INCISION/FLAP #2 OPEN/DEHISCED
VA-VAAES INCISION/FLAP #2 UNABLE
VA-VAAES INCISION/FLAP FLAP #2 IN PLACE
VA-VAAES INCISION/FLAP #2 WELL APPROX
VA-VAAES INCISION/FLAP #2 ASSESSED
VA-VAAES INCISION/FLAP #1 TEMP CHANGE
VA-VAAES INCISION/FLAP #1 MACERATION
VA-VAAES INCISION/FLAP #1 INDURATION
VA-VAAES INCISION/FLAP #1 ERYTHEMA
VA-VAAES INCISION/FLAP #1 INTACT
VA-VAAES INCISION/FLAP #1 OPEN/DEHISCED
VA-VAAES INCISION/FLAP #1 UNABLE
VA-VAAES INCISION/FLAP FLAP #1 IN PLACE
VA-VAAES INCISION/FLAP #1 WELL APPROX
VA-VAAES INCISION/FLAP #1 ASSESSED
VA-VAAES INCISION & FLAP ASSESS
SKIN INTEGRITY - WOUND
VA-VAAES WOUND TYPE PARTIAL THICKNESS
VA-VAAES WOUND TYPE FULL THICKNESS
SKIN INTEGRITY - WOUND ADDITIONAL
VAAES PI PERI WOUND ERYTHEMA NON-BLANCH
PI PERI WOUND INTACT
PI PERI WOUND USUAL FOR ETHNICITY
VA-VAAES OTHER WOUND #5 PERI WOUND - OTHER
VA-VAAES OTHER WOUND #5 PERI WOUND - SKIN TEMPERATURE
PI PERI WOUND FLUCTUANCE
PI PERI WOUND INDURATION
PI PERI WOUND CREPITUS
PI PERI WOUND MACERATION
PI PERI WOUND ERYTHEMA
VA-VAAES SKIN PI WOUND BED ESCHAR
VA-VAAES SKIN PI WOUND BED EPITH
VA-VAAES OTHER WOUND #5 WOUND BED TISSUE TYPE - OTHER
VA-VAAES SKIN PI WOUND BED NECROTIC
VA-VAAES SKIN PI WOUND BED SLOUGH
VA-VAAES SKIN PI WOUND BED GRANULATION
OTHER WOUND DRAINAGE VISCOUS

OTHER WOUND DRAINAGE SERO-PURULENT
OTHER WOUND DRAINAGE SEROUS
OTHER WOUND DRAINAGE SEROSANGUINEOUS
OTHER WOUND DRAINAGE PURULENT
OTHER WOUND DRAINAGE SANGUINEOUS
VA-VAAES OTHER WOUND DRAIN COLOR RED
VA-VAAES OTHER WOUND DRAIN COLOR AMBER
OTHER WOUND DRAINAGE COLOR GREEN
OTHER WOUND DRAINAGE COLOR TAN
OTHER WOUND DRAINAGE COLOR COLORLESS
OTHER WOUND DRAINAGE AMOUNT SCANT
OTHER WOUND DRAINAGE AMOUNT HEAVY
OTHER WOUND DRAINAGE AMOUNT MODERATE
OTHER WOUND DRAINAGE AMOUNT SMALL
OTHER WOUND WOUND DRAINAGE
OTHER WOUND DRAINAGE AMOUNT NONE
SKIN INTEGRITY - WOUND FIFTH
VA-VAAES OTHER WOUND #4 PERI WOUND - OTHER
VA-VAAES OTHER WOUND #4 PERI WOUND - SKIN TEMPERATURE
VA-VAAES OTHER WOUND #4 WOUND BED TISSUE TYPE - OTHER
SKIN INTEGRITY - WOUND FOURTH
VA-VAAES OTHER WOUND #3 PERI WOUND - OTHER
VA-VAAES OTHER WOUND #3 PERI WOUND - SKIN TEMPERATURE
VA-VAAES OTHER WOUND #3 WOUND BED TISSUE TYPE - OTHER
SKIN INTEGRITY - WOUND THIRD
VA-VAAES OTHER WOUND #2 PERI WOUND - OTHER
VA-VAAES OTHER WOUND #2 PERI WOUND - SKIN TEMPERATURE
VA-VAAES OTHER WOUND #2 WOUND BED TISSUE TYPE - OTHER
SKIN INTEGRITY - WOUND SECOND
VA-VAAES OTHER WOUND #1 PERI WOUND - OTHER
VA-VAAES OTHER WOUND #1 PERI WOUND - SKIN TEMPERATURE
VA-VAAES OTHER WOUND #1 WOUND BED TISSUE TYPE - OTHER
VA-VAAES SKIN INSPECTION PU/PI NEW
SKIN HIGH RISK FACTORS [C]
MED DEVICE PI MUCOSAL MEMBRANE FIFTH
PRESSURE ULCER STAGE [C]
MUCOSAL MEM PI #5
MED DEVICE PI MUCOSAL MEMBRANE FOURTH
MUCOSAL MEM PI #4
MED DEVICE PI MUCOSAL MEMBRANE THIRD
MUCOSAL MEM PI #3
MED DEVICE PI MUCOSAL MEMBRANE SECOND
MUCOSAL MEM PI #2
MED DEVICE PI MUCOSAL MEMBRANE

MUCOSAL MEM PI
MED DEVICE PI UNSTAGEABLE FIFTH
DRAINAGE APPEARANCE VISCOUS
DRAINAGE APPEARANCE SERO-PURULENT
DRAINAGE APPEARANCE PURULENT
DRAINAGE APPEARANCE SANGUINEOUS
DRAINAGE APPEARANCE SEROSANGUINEOUS
DRAINAGE APPEARANCE SEROUS
VA-VAAES DRAINAGE COLOR RED
VA-VAAES DRAINAGE COLOR AMBER
DRAINAGE COLOR GREEN
DRAINAGE COLOR TAN
DRAINAGE COLOR COLORLESS
DRAINAGE AMOUNT SCANT
DRAINAGE AMOUNT HEAVY
DRAINAGE AMOUNT MODERATE
DRAINAGE AMOUNT SMALL
WOUND DRAINAGE
DRAINAGE AMOUNT NONE
VA-VAAES PI UNSTAGEABLE #5 PERI WOUND - OTHER
VA-VAAES PI UNSTAGEABLE #5 PERI WOUND - SKIN TEMPERATURE
VA-VAAES PI UNSTAGEABLE #5 WOUND BED TISSUE TYPE - OTHER
UNABLE TO STAGE (PU) #5
MED DEVICE PI UNSTAGEABLE FOURTH
VA-VAAES PI UNSTAGEABLE #4 PERI WOUND - OTHER
VA-VAAES PI UNSTAGEABLE #4 PERI WOUND - SKIN TEMPERATURE
VA-VAAES PI UNSTAGEABLE #4 WOUND BED TISSUE TYPE - OTHER
UNABLE TO STAGE (PU) #4
MED DEVICE PI UNSTAGEABLE THIRD
VA-VAAES PI UNSTAGEABLE #3 PERI WOUND - OTHER
VA-VAAES PI UNSTAGEABLE #3 PERI WOUND - SKIN TEMPERATURE
VA-VAAES PI UNSTAGEABLE #3 WOUND BED TISSUE TYPE - OTHER
UNABLE TO STAGE (PU) #3
MED DEVICE PI UNSTAGEABLE SECOND
VA-VAAES PI UNSTAGEABLE #2 PERI WOUND - OTHER
VA-VAAES PI UNSTAGEABLE #2 PERI WOUND - SKIN TEMPERATURE
VA-VAAES PI UNSTAGEABLE #2 WOUND BED TISSUE TYPE - OTHER
UNABLE TO STAGE (PU) #2
MED DEVICE PI UNSTAGEABLE
VA-VAAES PI UNSTAGEABLE #1 PERI WOUND - OTHER
VA-VAAES PI UNSTAGEABLE #1 PERI WOUND - SKIN TEMPERATURE
VA-VAAES PI UNSTAGEABLE #1 WOUND BED TISSUE TYPE - OTHER
UNABLE TO STAGE (PU)
VA-VAAES PI STAGE 4 #1 PERI WOUND - OTHER

VA-VAAES PI STAGE 4 #1 PERI WOUND - SKIN TEMPERATURE
VA-VAAES PI STAGE 4 #1 WOUND BED TISSUE TYPE - OTHER
VA-VAAES PI STAGE 4 #5 PERI WOUND - OTHER
VA-VAAES PI STAGE 4 #5 PERI WOUND - SKIN TEMPERATURE
VA-VAAES PI STAGE 4 #5 WOUND BED TISSUE TYPE - OTHER
MED DEVICE PI STAGE 4 FIFTH
PRESENCE OF TENDON
PRESENCE OF CARTILAGE
PRESENCE OF BONE
VANOD UNDERMINING
VANOD TUNNELING
WOUND CLOCKFACE POSITION
STAGE IV #5
VA-VAAES PI STAGE 4 #4 PERI WOUND - OTHER
VA-VAAES PI STAGE 4 #4 PERI WOUND - SKIN TEMPERATURE
VA-VAAES PI STAGE 4 #4 WOUND BED TISSUE TYPE - OTHER
MED DEVICE PI STAGE 4 FOURTH
STAGE IV #4
VA-VAAES PI STAGE 4 #3 PERI WOUND - OTHER
VA-VAAES PI STAGE 4 #3 PERI WOUND - SKIN TEMPERATURE
VA-VAAES PI STAGE 4 #3 WOUND BED TISSUE TYPE - OTHER
MED DEVICE PI STAGE 4 THIRD
STAGE IV #3
VA-VAAES PI STAGE 4 #2 PERI WOUND - OTHER
VA-VAAES PI STAGE 4 #2 PERI WOUND - SKIN TEMPERATURE
VA-VAAES PI STAGE 4 #2 WOUND BED TISSUE TYPE - OTHER
MED DEVICE PI STAGE 4 SECOND
STAGE IV #2
MED DEVICE PI STAGE 4
STAGE IV
VA-VAAES PI STAGE 3 #1 PERI WOUND - OTHER
VA-VAAES PI STAGE 3 #1 PERI WOUND - SKIN TEMPERATURE
VA-VAAES PI STAGE 3 #1 WOUND BED TISSUE TYPE - OTHER
VA-VAAES PI STAGE 3 #5 PERI WOUND - OTHER
VA-VAAES PI STAGE 3 #5 PERI WOUND - SKIN TEMPERATURE
VA-VAAES PI STAGE 3 #5 WOUND BED TISSUE TYPE - OTHER
MED DEVICE PI STAGE 3 FIFTH
STAGE III #5
VA-VAAES PI STAGE 3 #4 PERI WOUND - OTHER
VA-VAAES PI STAGE 3 #4 PERI WOUND - SKIN TEMPERATURE
VA-VAAES PI STAGE 3 #4 WOUND BED TISSUE TYPE - OTHER
MED DEVICE PI STAGE 3 FOURTH
STAGE III #4
VA-VAAES PI STAGE 3 #3 PERI WOUND - OTHER

VA-VAAES PI STAGE 3 #3 PERI WOUND - SKIN TEMPERATURE
 VA-VAAES PI STAGE 3 #3 WOUND BED TISSUE TYPE - OTHER
 MED DEVICE PI STAGE 3 THIRD
 STAGE III #3
 VA-VAAES PI STAGE 3 #2 PERI WOUND - OTHER
 VA-VAAES PI STAGE 3 #2 PERI WOUND - SKIN TEMPERATURE
 VA-VAAES PI STAGE 3 #2 WOUND BED TISSUE TYPE - OTHER
 MED DEVICE PI STAGE 3 SECOND
 STAGE III #2
 MED DEVICE PI STAGE 3
 STAGE III
 VA-VAAES PI STAGE 2 #1 PERI WOUND - OTHER
 VA-VAAES PI STAGE 2 #1 PERI WOUND - SKIN TEMPERATURE
 VA-VAAES PI STAGE 2 #1 WOUND BED TISSUE TYPE -HYPER-PIGMENTATION
 VA-VAAES PI STAGE 2 #1 WOUND BED TISSUE TYPE - OTHER
 VA-VAAES PI STAGE 2 #1 WOUND BED TISSUE TYPE - RUPTURED BLISTER
 VA-VAAES PI STAGE 2 #1 WOUND BED TISSUE TYPE - INTACT BLISTER
 VA-VAAES PI STAGE 2 #1 WOUND BED TISSUE TYPE - RED
 VA-VAAES PI STAGE 2 #1 WOUND BED TISSUE TYPE - PINK
 VA-VAAES PI STAGE 2 #1 WOUND BED TISSUE TYPE - HYPO-PIGMENTATION
 VA-VAAES PI STAGE 2 #5 PERI WOUND - OTHER
 VA-VAAES PI STAGE 2 #5 PERI WOUND - SKIN TEMPERATURE
 VA-VAAES PI STAGE 2 #5 WOUND BED TISSUE TYPE -HYPER-PIGMENTATION
 VA-VAAES PI STAGE 2 #5 WOUND BED TISSUE TYPE - OTHER
 VA-VAAES PI STAGE 2 #5 WOUND BED TISSUE TYPE - RUPTURED BLISTER
 VA-VAAES PI STAGE 2 #5 WOUND BED TISSUE TYPE - INTACT BLISTER
 VA-VAAES PI STAGE 2 #5 WOUND BED TISSUE TYPE - RED
 VA-VAAES PI STAGE 2 #5 WOUND BED TISSUE TYPE - PINK
 VA-VAAES PI STAGE 2 #5 WOUND BED TISSUE TYPE - HYPO-PIGMENTATION
 MED DEVICE PI STAGE 2 FIFTH
 STAGE II #5
 VA-VAAES PI STAGE 2 #4 PERI WOUND - OTHER
 VA-VAAES PI STAGE 2 #4 PERI WOUND - SKIN TEMPERATURE
 VA-VAAES PI STAGE 2 #4 WOUND BED TISSUE TYPE-HYPER-PIGMENTATION
 VA-VAAES PI STAGE 2 #4 WOUND BED TISSUE TYPE - OTHER
 VA-VAAES PI STAGE 2 #4 WOUND BED TISSUE TYPE - RUPTURED BLISTER
 VA-VAAES PI STAGE 2 #4 WOUND BED TISSUE TYPE - INTACT BLISTER
 VA-VAAES PI STAGE 2 #4 WOUND BED TISSUE TYPE - RED
 VA-VAAES PI STAGE 2 #4 WOUND BED TISSUE TYPE - PINK
 VA-VAAES PI STAGE 2 #4 WOUND BED TISSUE TYPE - HYPO-PIGMENTATION
 MED DEVICE PI STAGE 2 FOURTH
 STAGE II #4
 VA-VAAES PI STAGE 2 #3 PERI WOUND - OTHER
 VA-VAAES PI STAGE 2 #3 PERI WOUND - SKIN

TEMPERATURE
VA-VAAES PI STAGE 2 #3 WOUND BED TISSUE TYPE
-HYPER-PIGMENTATION
VA-VAAES PI STAGE 2 #3 WOUND BED TISSUE TYPE -
OTHER
VA-VAAES PI STAGE 2 #3 WOUND BED TISSUE TYPE -
RUPTURED BLISTER
VA-VAAES PI STAGE 2 #3 WOUND BED TISSUE TYPE -
INTACT BLISTER
VA-VAAES PI STAGE 2 #3 WOUND BED TISSUE TYPE - RED
VA-VAAES PI STAGE 2 #3 WOUND BED TISSUE TYPE - PINK
VA-VAAES PI STAGE 2 #3 WOUND BED TISSUE TYPE - HYPO-PIGMENTATION
MED DEVICE PI STAGE 2 THIRD
STAGE II #3
VA-VAAES PI STAGE 2 #2 PERI WOUND - OTHER
VA-VAAES PI STAGE 2 #2 PERI WOUND - SKIN TEMPERATURE
VA-VAAES PI STAGE 2 #2 WOUND BED TISSUE TYPE-HYPER-PIGMENTATION
VA-VAAES PI STAGE 2 #2 WOUND BED TISSUE TYPE - OTHER
VA-VAAES PI STAGE 2 #2 WOUND BED TISSUE TYPE - RUPTURED BLISTER
VA-VAAES PI STAGE 2 #2 WOUND BED TISSUE TYPE - INTACT BLISTER
VA-VAAES PI STAGE 2 #2 WOUND BED TISSUE TYPE - RED
VA-VAAES PI STAGE 2 #2 WOUND BED TISSUE TYPE - PINK
VA-VAAES PI STAGE 2 #2 WOUND BED TISSUE TYPE - HYPO-PIGMENTATION
MED DEVICE PI STAGE 2 SECOND STAGE II #2
MED DEVICE PI STAGE 2 STAGE II
MED DEVICE PI STAGE 1 FIFTH STAGE I #5
MED DEVICE PI STAGE 1 FOURTH STAGE I #4
MED DEVICE PI STAGE 1 THIRD STAGE I #3
MED DEVICE PI STAGE 1 SECOND STAGE I #2
MED DEVICE PI STAGE 1 STAGE I
MED DEVICE PI DEEP TISSUE FIFTH
VA-VAAES PI DEEP TISSUE #5 PERI WOUND - OTHER
VA-VAAES PI DEEP TISSUE #5 PERI WOUND - SKIN TEMPERATURE
VA-VAAES PI DEEP TISSUE #5 WOUND BED TISSUE TYPE - OTHER
SUSPECTED DEEP TISSUE INJURY #5
MED DEVICE PI DEEP TISSUE FOURTH
VA-VAAES PI DEEP TISSUE #4 PERI WOUND - OTHER
VA-VAAES PI DEEP TISSUE #4 PERI WOUND - SKIN TEMPERATURE
VA-VAAES PI DEEP TISSUE #4 WOUND BED TISSUE TYPE - OTHER
SUSPECTED DEEP TISSUE INJURY #4
MED DEVICE PI DEEP TISSUE THIRD
VA-VAAES PI DEEP TISSUE #3 PERI WOUND - OTHER
VA-VAAES PI DEEP TISSUE #3 PERI WOUND - SKIN TEMPERATURE
VA-VAAES PI DEEP TISSUE #3 WOUND BED TISSUE TYPE - OTHER

SUSPECTED DEEP TISSUE INJURY #3
MED DEVICE PI DEEP TISSUE SECOND
VA-VAAES PI DEEP TISSUE #2 PERI WOUND - OTHER
VA-VAAES PI DEEP TISSUE #2 PERI WOUND - SKIN TEMPERATURE
VA-VAAES PI DEEP TISSUE #2 WOUND BED TISSUE TYPE - OTHER
SUSPECTED DEEP TISSUE INJURY #2
MED DEVICE PI DEEP TISSUE
VA-VAAES PI DEEP TISSUE #1 PERI WOUND - OTHER
VA-VAAES PI DEEP TISSUE #1 PERI WOUND - SKIN TEMPERATURE
VA-VAAES PI DEEP TISSUE #1 WOUND BED TISSUE TYPE - OTHER
SUSPECTED DEEP TISSUE INJURY
PRESSURE ULCER
SKIN INTEGRITY - FISSURE
SKIN PROBLEM - TEAR
SKIN PROBLEM - SCRATCHES
VA-VAAES SKIN INTEGRITY LACERATION
SKIN INTEGRITY - HEMATOMA
SKIN PROBLEM - ERYTHEMIA NONBLANCHABLE
SKIN PROBLEM - CHAFING
SKIN PROBLEM - OTHER
VA-VAAES SKIN INTEGRITY RASH OTHER
VA-VAAES SKIN INTEGRITY RASH VESICLES
VA-VAAES SKIN INTEGRITY RASH LOCALIZED
VA-VAAES SKIN INTEGRITY RASH HIVES
VA-VAAES SKIN INTEGRITY RASH GENERALIZED
VA-VAAES SKIN INTEGRITY RASH DIFFUSE
SKIN PROBLEM - RASH
SKIN INTEGRITY - MASD
SKIN INTEGRITY - MACERATION
SKIN INTEGRITY - EXCORIATION
SKIN INTEGRITY - ERYTHEMA
SKIN INTEGRITY - DENUDED
SKIN PROBLEMS - BURN
SKIN PROBLEM - BRUISING
VA-VAAES SKIN INTEGRITY BLISTER OTHER
VA-VAAES SKIN INTEGRITY BLISTER WEEPING
VA-VAAES SKIN INTEGRITY BLISTER UNROOFED
VA-VAAES SKIN INTEGRITY BLISTER INTACT
SKIN INTEGRITY - BLISTER
VA-VAAES SKIN INTEGRITY BITE TICK
VA-VAAES SKIN INTEGRITY BITE OTHER
VA-VAAES SKIN INTEGRITY BITE SPIDER
VA-VAAES SKIN INTEGRITY BITE SNAKE
VA-VAAES SKIN INTEGRITY BITE INSECT

VA-VAAES SKIN INTEGRITY BITE HUMAN
VA-VAAES SKIN INTEGRITY BITE FELINE
VA-VAAES SKIN INTEGRITY BITE CANINE
SKIN INTEGRITY - BITE WOUND
VA-VAAES SKIN INTEGRITY ABRASION
VA-VAAES SKIN INTEGRITY INTACT
PRESSURE ULCER REASSESSMENT [C]
NEW PRESSURE ULCER
PRESSURE ULCER REASSESS - CHANGE
SKIN HIGH RISK OTHER
POTENTIAL COMPROMISED NUTRITION
NEUROPATHY
MEDICAL DEVICE
MALNUTRITION DIAGNOSIS
VASCULAR SURGERY OR VASCULAR DISEASE
VA-VAAES PI RISK PREVIOUS HISTORY
AMPUTEE
AGE OVER 75
VA-VAAES SKIN RISK FACTOR NO
PRESSURE ULCER-EDUCATION [C]
EDUCATION MATERIALS ON ULCER PREVENTION
EDUCATE IMPORTANCE OF CHANGING POSITION
PROVIDE EDUCATION REGARDING TX PLAN
PROVIDE EDUCATION ON CAUSE/PREVENTION
PU EDUCATION-OTHER
VA-VAAES PRESSURE INJURY INTERVENTIONS [C]
SKIN - EDUCATION
PRESSURE ULCER-FRICTION/SHEAR [C]
PU FRICTION/SHEAR-OTHER
USE BED TRAPEZE OR PULL SHEET
WHEN HOB ELEVATED RAISE KNEE
HOB BELOW 30 DEGREES WHEN NOT EATING
ELEVATE HOB FOR MEALS
SKIN REDUCE FRICTION AND SHEAR
PRESSURE ULCER-NUTRITION [C]
PU NUTRITION-OTHER
TUBE FEEDINGS AS ORDERED
TRAY SET UP AND ASSISTANCE
PROVIDE/ENCOURAGE ORAL CARE AS NEEDED
OFFER ORDERED SUPPLEMENTS
OFFER LIQUIDS Q2H WHEN TURNING
MONITOR FLUID/FOOD INTAKE
ENCOURAGE EATING AND ASSIST WITH MEALS
SKIN MANAGE NUTRITION

PRESSURE ULCER-MOISTURE [C]
SCHEDULED TOILETING
PU MOISTURE-OTHER
OFFER BEDPAN/URINAL
NO MORE THAN 1 LINEN LAYER
MAINTAIN CLEAN DRY SKIN
INSTRUCT PT/FAMILY TO REQUEST ASSISTANCE
PROTECTIVE BARRIER OINTMENT
FECAL COLLECTOR
CONDOM CATHETER
SKIN MANAGE MOISTURE
PRESSURE ULCER-REMOBILIZE [C]
ROM EXERCISES
PU REMOBILIZE-OTHER
LIMIT SITTING OOB TO 2 HR PERIODS
REPOSITION IN CHAIR
ENCOURAGE ACTIVITY AS TOLERATED
SKIN - MAXIMAL REMOBILIZATION
PRESSURE ULCER-PRESSURE REDUCING [C]
SPECIALTY BED/SURFACE
PU PRESSURE REDUCING-OTHER
WHEELCHAIR CUSHION
TURN AND REPOSITION Q2H
FREQUENT POSITION CHANGES
ELEVATE HEELS
TURN TO SIDE LESS THAN 30 DEGREES
HEEL/ELBOW PADS
SKIN - PRESSURE-REDISTRIBUTION MEASURES
VA-VAAES PRESSURE INJURY INTERVENTION
VA-VAAES PRESSURE INJURY INT NOT NEEDED
INDIVIDUALIZED REPOSITIONING
ENCOURAGE MEALS AND ASSIST AS NEEDED
VA-VAAES PRESSURE INJURY INT NO CHANGE
VA-VAAES SKIN INSPECTION/ASSESSMENT V4.1

REMINDER SPONSOR

OFFICE OF NURSING SERVICE
MENTAL HEALTH SERVICES

REMINDER TERM

VA-REMINDER UPDATE_2_0_460
VA-VANOD SKIN INSPECTION ONLY
VA-BL NA USER CLASS
VA-BL LPN USER CLASS

HEALTH SUMMARY TYPE

VAAES PAST WOUNDS
VA-PRESSURE ULCER
VA-PU INTERVENTIONS
VA-VANOD BRADEN SCALE

HEALTH SUMMARY OBJECTS

VAAES PAST WOUNDS
VA-PRESSURE ULCER (TIU)
VA-PU INTERVENTIONS (TIU)
VA-VANOD BRADEN SCALE

TIU DOCUMENT DEFINITION

VAAES PAST WOUNDS
PRESSURE ULCER
PU INTERVENTIONS
VA-VANOD BRADEN SCALE

REMINDER DIALOG

VAAES SKIN INSPECTION/ASSESSMENT

Pre-Installation

Using the VistA menu, Reminder Dialog Manager, Reminder Dialogs, change view to Element View. Inquire to the following elements to check for local quick orders in the finding item field:

VAL-VANOD OI WOUND
VAL-VANOD OI NUTRITION CONSULT

Install Details

This update is being distributed as a web host file. The address for the host file is:

https://REDACTED/UPDATE_2_0_460.PRD

The file will be installed using Reminder Exchange, programmer access is not required.

Installation:

=====

This update can be loaded with users on the system. Installation will take less than 30 minutes.

Install Example

To Load the Web Host File. Navigate to Reminder exchange in Vista


```

+      + Next Screen  - Prev Screen  ?? More Actions
CFE Create Exchange File Entry      LHF Load Host File
CHF Create Host File                LMM Load MailMan Message
CMM Create MailMan Message          LR  List Reminder Definitios
DFE Delete Exchange File Entry      LWH Load Web Host File
IFE Install Exchange File Entry      RI  Reminder Definition Inqy
IH  Installation History             RP  Repack
Select Action: Next Screen// LWH  Load Web Host File

```

At the **Select Action:** prompt, enter **LWH** for Load Web Host File

At the **Input the url for the .prd file:** prompt, type the following web address:

https://REDACTED/UPDATE_2_0_460.PRD

You should see a message at the top of your screen that the file successfully loaded.

Search for and locate the reminder exchange entry titled UPDATE_2_0_460 VA-VAAES SKIN INSPECTION/ASSESSMENT UPDATE.

```

+Item  Entry                               Source                               Date Packed
 743  UPDATE_2_0_460 VA-VAAES SKIN          [REDACTED]                          [REDACTED]
      INSPECTION/ASSESSMENT UPDATE
 744  UPDATE_2_0_461 VA-EXTENDED          [REDACTED]                          [REDACTED]
      REALITY REHABILITATION
      TEMPLATE
 745  UPDATE_2_0_462 VA-ONCOLOGY          [REDACTED]                          [REDACTED]
      COLON CANCER CLINICAL PATHWAYS
 746  UPDATE_2_0_47 VA-HEPATITIS B       [REDACTED]                          [REDACTED]
      IMMUNIZATIONS
 747  UPDATE_2_0_48 MISC UPDATE 6        [REDACTED]                          [REDACTED]
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE Create Exchange File Entry      LHF Load Host File
CHF Create Host File                LMM Load MailMan Message
CMM Create MailMan Message          LR  List Reminder Definitions
DFE Delete Exchange File Entry      LWH Load Web Host File
IFE Install Exchange File Entry      RI  Reminder Definition Inquiry
IH  Installation History             RP  Repack
Select Action: Next Screen// IFE  Install Exchange File Entry
Enter a list or range of numbers (1-970): 743

```

At the **Select Action** prompt, enter **IFE** for Install Exchange File Entry

Enter the number that corresponds with your entry titled **UPDATE_2_0_460 VA-VAAES SKIN INSPECTION/ASSESSMENT UPDATE** (in this example it is entry 743, it will vary by site). The exchange file date should be 07/25/2024.

```
Package Version: 2.0P65

Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
  VAAES SKIN INSPECTION/ASSESSMENT

REMINDER TERM
  VA-REMINDER UPDATE_2_0_460

Keywords:

Components:

+      Enter ?? for more actions                                >>>
IA   Install all Components                                IS   Install Selected Component
Select Action: Next Screen// IA
```

At the **Select Action** prompt, type **IA** for Install all Components and hit enter.

Select Action: Next Screen// **IA Install all Components**

You will see several prompts:

Choose **I to Install** all new components

If you are prompted to overwrite any existing components, choose **S to Skip**.

You will be prompted to install the dialog component:

```
Packed reminder dialog: VAAES SKIN INSPECTION/ASSESSMENT
```

Item	Seq.	Dialog Findings	Type	Exists
1		VAAES SKIN INSPECTION/ASSESSMENT	dialog	X
2	1	VA-VANOD VERSION NUMBER Finding: VA-VAAES SKIN INSPECTION/ASSESSMENT V3.2 (HEALTH FACTOR)	element	X
3	2	VA-VANOD MONITORING AND RISK ASSESSMENT Finding: *NONE*	group	X
4	2.5	VA-VANOD MONITORING AND RISK LINKS Finding: *NONE*	element	X
5	3	VA-VANOD GP INITIAL OR REASSESSMENT-G* Finding: *NONE*	group	
6	3.10	VA-VANOD GP SKIN ASSESSMENT (COMPLETE) Finding: *NONE*	group	X

```
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// IA
```

At the **Select Action** prompt, type **IA** to install the dialog – **VAAES SKIN INSPECTION/ASSESSMENT**.

Select Action: Next Screen// **IA Install All**

Install reminder dialog and all components with no further changes: **Y// ES**

```
Packed reminder dialog: VAAES SKIN INSPECTION/ASSESSMENT
VAAES SKIN INSPECTION/ASSESSMENT (reminder dialog) installed from exchange file.
```

Item	Seq.	Dialog Findings	Type	Exists
1		VAAES SKIN INSPECTION/ASSESSMENT	dialog	X
2	1	VA-VANOD VERSION NUMBER Finding: VA-VAAES SKIN INSPECTION/ASSESSMENT V3.2 (HEALTH FACTOR)	element	X
3	2	VA-VANOD MONITORING AND RISK ASSESSMENT Finding: *NONE*	group	X
4	2.5	VA-VANOD MONITORING AND RISK LINKS Finding: *NONE*	element	X
5	3	VA-VANOD GP INITIAL OR REASSESSMENT-G* Finding: *NONE*	group	X
6	3.10	VA-VANOD GP SKIN ASSESSMENT (COMPLETE) Finding: *NONE*	group	X

```
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// Q
```

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

```
Package Version: 2.0P65

Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
  VAAES SKIN INSPECTION/ASSESSMENT

REMINDER TERM
  VA-REMINDER UPDATE_2_0_460

Keywords:

Components:

+      + Next Screen  - Prev Screen  ?? More Actions  >>>
IA  Install all Components          IS  Install Selected Component
Select Action: Next Screen// Q
```

You will be returned to this screen. Type Q to quit.

Install Completed

Post Installation

1. **If you were not prompted during the install to replace with your local consult order, add your local quick orders as finding items to the elements listed below.**
Using the Vista menu, reminder dialog manager, reminder dialogs, change view to Element view and search for the following elements:
 - VAL-VANOD OI WOUND
 1. At the finding item prompt, add local quick order for consult for Wound in the finding item field.
 - VAL-VANOD OI NUTRITION CONSULT
 1. At the finding item prompt, add local quick order for consult order for nutrition in the finding item field.

2. **Determine if Pressure Ulcer Injury documentation display should remain enabled.** Per field request, the group VAL-GP VANOD PRIOR PRESSURE ULCER/INJURY was made a local group which will allow sites to disable this information from displaying in the template and printing to the final note if selected by clinical staff.

3. In CPRS, open a new VAAES SKIN INSPECTION/ASSESSMENT note and confirm the version number now displays as 4.1.

Install complete.